

Cynomolgus Monkey Serum IgG Depleted

Catalog No.	CS209A	Quantity:	50 ml
	CS209B		100 ml

Description: Prepared from mixed gender, pooled, frozen serum from *Macaca fascicularis*, using immobilized Protein A.

Source: *Macaca fascicularis*

Formulation: Frozen liquid

Residual IgG: <0.5 mg/ml Cynomolgus Monkey IgG

Storage & Stability: When stored at -80°C, product is stable for 3 years from date of delivery. Upon initial thaw, it is recommended to prepare working aliquots and store at -80°C. **Avoid repeated freeze-thaw cycles.**

Infectious Disease Testing: Non-human primate products are sourced from animals screened and found negative for Herpes-B Virus, SRV, SIV, STLV-1, Measles and TB

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

