

Cynomolgus Monkey Serum, IgG Depleted

Catalog No.	CS209A CS209B	Quantity:	50 ml 100 ml
Description:	Prepared from frozen cynomolgus monkey serum using immobilized Protein A.		
Source:	<i>Macaca fascicularis</i> serum		
Formulation:	Frozen liquid		
Residual IgG:	< 500 ng/ml Cynomolgus Monkey IgG, lot specific, typically < 100 ng/mg		
Storage & Stability:	Store at -80°C for up to 1 year. Upon initial thaw, prepare working aliquots and store at -80°C. Avoid repeated freeze-thaw cycles.		
Infectious Disease Testing:	Non-human primate products are sourced from animals screened and found negative for Herpes-B Virus, SRV, SIV, STLV-1, Measles and TB□		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



Cell Sciences®
65 Parker Street
Unit 11
Newburyport, MA 01950

Toll Free: 888-769-1246
Phone: 978-572-1070
Fax: 978-992-0298

E-mail: info@cellsciences.com
Website: www.cellsciences.com