

Cyno Monkey IgG Depleted Serum

Catalog No.	CS209A CS209B	Quantity:	50 ml 100 ml
Description:	Prepared from frozen cyno monkey serum using immobilized Protein A.		
Concentration:	<0.5 mg/ml Cyno Monkey IgG		
Purity:	Serum		
Formulation:	Liquid		
Storage & Stability:	When stored at -80°C, product is stable for 3 years from date of delivery. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

