

PLG

Native Equine Plasminogen

Catalog No.	CS213A	Quantity:	1.0 mg
	CS213B		10 mg

Description: Prepared from fresh horse plasma by immobilized lysine chromatography. Please inquire about the availability of active horse plasmin.

UniProt ID: P80010

Gene ID: 100056744

Source: Horse plasma

Molecular Weight: 92 kDa

Formulation: Liquid in 0.1M HEPES, 0.1M NaCl, pH 7.4

Purity: >95% by SDS-PAGE analysis

Concentration: ≥ 1.0 mg/ml, lot specific

Extinction Coefficient: $E^{0.1\%}_{280\text{nm}} = 1.69$

Storage & Stability: Store at -80°C for up to 1 year. Upon initial thaw, prepare working aliquots and store at -80°C . **Avoid repeated freeze-thaw cycles.**

1 Horse Plasminogen (3 μg) NR



NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

cellsciences.com



Cell Sciences®
65 Parker Street
Unit 11
Newburyport, MA 01950

Toll Free: 888-769-1246
Phone: 978-572-1070
Fax: 978-992-0298

E-mail: info@cellsciences.com
Website: www.cellsciences.com