

## Serpine 1

### Mouse PAI-1 (Biotin labeled active fraction)

<b>Catalog No.</b>	CS223A CS223B	<b>Quantity:</b>	0.5 mg 1.0 mg
<b>Alternate Names:</b>	PAI1, PAI-1, Planh1		
<b>Description:</b>	This active fraction of recombinant wild type mouse PAI-1 is biotin labeled on lysine residues. Typical preparations are ~50% active.		
<b>Concentration:</b>	2.5 mg/ml		
<b>Gene ID:</b>	18787		
<b>Source:</b>	<i>E. coli</i>		
<b>Molecular Weight:</b>	43 kDa		
<b>Formulation:</b>	Liquid in 0.05M Sodium Acetate; 0.1M NaCl; 1mM EDTA; pH 5.0		
<b>Purity:</b>	>95% by SDS-PAGE analysis		
<b>Storage &amp; Stability:</b>	When stored at -80°C, product is stable for 3 years from date of delivery. <b>Avoid repeated freeze-thaw cycles.</b>		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

