

## Native Rabbit IgG FITC pAb

<b>Catalog No.</b>	CS237A CS237B	<b>Quantity:</b>	1.0 mg 10 mg
<b>Description:</b>	Rabbit IgG purified from normal serum by immobilized Protein A then FITC labeled.		
<b>Concentration:</b>	10 mg/ml		
<b>Molecular Weight:</b>	160 kDa		
<b>Purity:</b>	>95% by SDS-PAGE analysis and preservative free		
<b>Specificity:</b>	This antibody has been selected as an isotype control as it has no known specificity.		
<b>Host:</b>	Rabbit		
<b>Source:</b>	Rabbit plasma		
<b>Isotype:</b>	IgG		
<b>Clone:</b>	Polyclonal		
<b>Conjugate:</b>	FITC		
<b>Formulation:</b>	15 mM PBS, pH 7.4, without preservative		
<b>Purification:</b>	Protein A affinity chromatography		
<b>Extinction Coefficient:</b>	$E^{0.1\%}_{280\text{nm}} = 1.36$		
<b>Applications:</b>	ELISA Western Blot		
<b>Storage &amp; Stability:</b>	Store at -80°C for up to 1 year. Upon initial thaw, prepare working aliquots and store at -80°C. <b>Avoid repeated freeze-thaw cycles.</b>		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

