

Mouse Anti-Human Factor V Clone 14G5 Biotin mAb

Catalog No.	CS372A CS372B	Quantity:	100 µg 1.0 mg
Description:	Biotin labeled mouse monoclonal antibody to human Factor V. IgG fraction purified by immobilized Protein G. Clone 14G5-G11, isotype IgG1Kappa		
Concentration:	1.0 mg/ml		
Gene ID:	2153		
Purity:	IgG fraction		
Clone:	14G5-G11		
Host:	Mouse		
Immunogen:	Full length native protein (purified) (Human).		
Isotype:	IgG1 kappa		
Molecular Weight:	160 kDa		
Source:	Hybridoma cell culture		
Formulation:	Liquid in 0.0 5M Sodium Phosphate + 0.1M NaC + 1 mM EDTA, pH 6.6		
Purification:	IgG fraction.		
Cross-Reactivity:	Human		
Applications:	ELISA, WB		
Storage & Stability:	When stored at -80°C, product is stable for 5 years from date of delivery. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

