

Mouse Cerebral Cortex Tissue Lysate

Catalog No.	CS378A	Quantity:	200 µg
	CS378B		400 µg

Description: A freshly-harvested Balb C mouse brain is dissected and the cerebral cortex is flash-frozen in liquid nitrogen and stored frozen until lysis. The tissue is ground and thoroughly homogenized at 4C with a freshly prepared modified Radio Immuno Precipitation Assay (RIPA) buffer containing strong protease and phosphatase inhibitor cocktails. The final homogenate is incubated with DNase then clarified by centrifugation at 10,000 g at 4C and stored at -80C. Each lot is quality controlled by determining the protein concentration using a BCA assay and analyzed with SDS-PAGE using Coomassie blue for visualization. Sharp, clear protein bands indicate a quality protein lysate.

Concentration: 2.5 mg/ml

Source: Mouse brain

Formulation: Liquid in 0.1 M Tris-HCl + 0.15 M NaCl + 0.1% SDS + 0.5% Na Deoxycholate + 1% Triton X-100, pH 7.4

Storage & Stability: When stored at -80°C, product is stable for 5 years from date of delivery. **Avoid repeated freeze-thaw cycles.**



1. Mouse Cerebral Cortex Lysate (10 µg) R

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



Cell Sciences®
480 Neponset Street
Bldg 12A
Canton, MA 02021

Toll Free: 888-769-1246
Phone: 781-828-0610
Fax: 781-828-0542

E-mail: info@cellsciences.com
Website: www.cellsciences.com