

KLK3

Sheep Anti-Human Prostate Specific Antigen Antiserum pAb

Catalog No.	CS398A CS398B	Quantity:	1.0 ml 10 ml
Description:	Polyclonal antiserum (host rabbit) to human PSA.		
Concentration:	Whole Antiserum		
Gene ID:	354		
Purity:	Serum		
Host:	Sheep		
Immunogen:	Full length native protein (purified) (Human).		
Isotype:	IgG		
Formulation:	Frozen liquid		
Cross-Reactivity:	Human		
Applications:	ELISA, WB		
Storage & Stability:	When stored at -80°C, product is stable for 5 years from date of delivery. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

