

Sheep Anti-Swine IgG Protein G Purified pAb

Catalog No.	CS405A CS405B	Quantity:	1 g 2 g
Description:	Sheep Anti-Swine IgG Protein G Purified polyclonal antibody.		
Concentration:	10 mg/ml		
Purity:	>95 percent pure by SDS-PAGE		
Molecular Weight:	160 kDa		
Host:	Sheep		
Isotype:	IgG		
Target:	This antibody has been selected as an isotype control as it has no known specificity.		
Formulation:	Frozen liquid in 0.02 M Sodium Phosphate + 0.15 M NaCl, pH 7.4. Preservative free.		
Purification:	Purified from normal serum by immobilized Protein G.		
Applications:	ELISA WB		
Storage & Stability:	When stored at -80°C, product is stable for 5 years from date of delivery. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

