

ADGRB3

Human Adhesion G Protein-Coupled Receptor B3

Catalog No.	CSH3048MP	Quantity:	10 mg
	CSH3048PR		50 µg

Alternate Names: Brain-Specific Angiogenesis Inhibitor 1, BAI1, GDAIF

Description: Adhesion G protein-coupled receptors (adhesion GPCRs) are a class of 33 human protein receptors with a broad distribution in embryonic and larval cells, cells of the reproductive tract, neurons, leukocytes, and a variety of tumours. The defining feature of adhesion GPCRs that distinguishes them from other GPCRs is their hybrid molecular structure. The extracellular region of adhesion GPCRs can be exceptionally long and contain a variety of structural domains that are known for the ability to facilitate cell and matrix interactions. ADGRB3 is postulated to be a member of the secretin receptor family and acts as an angiogenesis inhibitor.

The receptor is available in the following formats: stable over-expression cell line, membrane preparation, or purified receptor in HEK293 or CHO. Various tagged versions are available.

Gene ID: 577

UniProtKB: O60242

Format: Cell line, membrane preparation, or purified protein

Source: HEK 293 or CHO cells

Characterization: Expression of receptor was verified by immunostaining. Receptor demonstrates biological activity when tested in a radioligand assay.

Affinity Tag Options: Receptor construct: A3 is 2X Twin-Strep tagged and HIS tagged

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.