

## Native Bovine Gamma Globulin

<b>Catalog No.</b>	CSI10305	<b>Quantity:</b>	10 g
<b>Description:</b>	Gamma globulins are a class of globulins, identified by their position after serum protein electrophoresis. The most significant gamma globulins are immunoglobulins ("Igs"), more commonly known as antibodies, although some Igs are not gamma globulins, and some gamma globulins are not Igs.		
<b>Source:</b>	Bovine serum		
<b>Formulation:</b>	Lyophilized in 0.01M PBS buffer, pH 7.2.		
<b>Purity:</b>	>99% pure		
<b>Reconstitution:</b>	Reconstitute in laboratory grade water, or buffer of your choice.		
<b>Applications:</b>	User optimized		
<b>Storage &amp; Stability:</b>	Store at 4° C.		

**NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.**