

Mouse Immunoglobulin G2b kappa Isotype Control Clone MG2b-57 LE/NA Antibody

Catalog No.	CSI11943 CSI11944	Quantity:	50 µg 500 µg
Description:	The isotype of MG2b-57 immunoglobulin is mouse IgG2b,κ. This antibody was chosen as an isotype control after screening on a variety of resting, activated, live, and fixed mouse, rat and human tissues.		
Concentration:	1.0 mg/ml		
Purity:	>95% by SDS-PAGE analysis.		
Immunogen:	Trinitrophenol + KLH		
Isotype:	Mouse IgG2b, κ		
Clone:	MG2b-57		
Formulation:	Phosphate-buffered solution, pH 7.2, containing no preservative. 0.2 µm filter sterilized. Endotoxin level is < 0.1 EU/µg of the protein (< 0.01 ng/µg of the protein) as determined by the LAL test.		
Purification:	The (Low Endotoxin, Azide-Free) antibody was Purified by affinity chromatography.		
Applications:	FC, ICFC, WB, IP, ICC, IF, FA		
Recommended Usage:	Each lot of this antibody is quality control tested by immunofluorescent staining with flow cytometric analysis as negative control, and the purity is greater than 95% by SDS-PAGE. Use at concentrations comparable to those of the specific antibody of interest.		
Storage & Stability:	The antibody solution should be stored undiluted at 4 °C. This LE/NA solution contains no preservative; handle under aseptic conditions		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

