

IgG

Native Mouse Immunoglobulin G

Catalog No.	CSI11978	Quantity:	10 mg
Alternate Names:	Mouse IgG Protein, Immunoglobulin G, Mouse IgG, Mouse IgG, Purified IgG		
Description:	Purified Mouse IgG - Immunoglobulin G (IgG) are antibody molecules. Each IgG is composed of four peptide chains - two heavy chains and two light chains. Each IgG has two antigen binding sites. Other Immunoglobulins may be described in terms of polymers with the IgG structure considered the monomer.		
Formulation:	Lyophilized from 10 mM PO4 buffer with 0.15 M NaCl, pH 7.2 and no preservatives.		
Purity:	>97% pure		
Reconstitution:	Reconstitute with laboratory grade water.		
Applications:	User optimized		
Storage & Stability:	Aliquot and store at -20°C. Avoid repeated Freeze/Thaw cycles		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.