

Human Immunoglobulin G (IgG)

Catalog No.	CSI11983A	Quantity:	1.0 gram
	CSI11983B		5.0 gram

Description: Human Immunoglobulin G (IgG) is a monomeric immunoglobulin, built of two heavy chains and two light chains. Each molecule has two antigen binding sites. Human Immunoglobulin G (IgG) is the most abundant immunoglobulin and is approximately equally distributed in blood and in tissue liquids.

Source: Human Plasma

Formulation: Lyophilized

Purity: >98% by SDS PAGE

Storage & Stability: Store at 2-4°C.

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

