

Recombinant Treponema pallidum p15 Protein

Catalog No.	CSI13230	Quantity:	0.5 mg
Alternate Names:	T.pallidum p15 protein		
Description:	Purified recombinant Treponema pallidum p15 protein - Treponema pallidum, the causative spirochete bacteria in syphilis, is a highly invasive organism that is capable of infecting all body tissues. The bacterial cells have a coat of glycosaminoglycans and the interior, unusually, has a very high lipid level including cardiolipin.		
Source:	<i>E. coli</i>		
Species:	Viral		
Formulation:	0.1M PBS, pH 7.4 with 0.01% NaN ₃ . Precaution: Sodium azide is a poisonous and hazardous substance which should be handled by trained staff only.		
Purity:	>99% pure by SDS PAGE		
Applications:	User optimized		
Storage & Stability:	Aliquot and store at -70°C or lower. Avoid repeated freeze/thaw cycles		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.