

ALB

BSA, Rabbit Anti-Bovine Serum Albumin pAb

Catalog No.	CSI19719A	Quantity:	1 mL
	CSI19719B		10 mL

Description: Bovine Serum Albumin (BSA), also known as "Fraction V", is a serum albumin protein that has numerous biochemical applications including ELISAs, immunoblots, and immunohistochemistry. It is also used as a nutrient in cell and microbial culture. BSA is also commonly used to determine the quantity of other proteins, by comparing an unknown quantity of protein to known amounts of BSA. BSA is used because of its stability and its lack of effect in many biochemical reactions.

Concentration: Whole Antiserum

Volume: 1.0 mL

Host: Rabbit Serum

Formulation: Frozen Liquid

Applications: ELISA, WB.
The optimal concentration should be determined by the user for each specific application.

Storage & Stability: Store at -70°C.

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

