

## Serpine1 Mouse PAI-1 Depleted Plasma (Na Citrate)

<b>Catalog No.</b>	CSI19894A CSI19894B	<b>Quantity:</b>	10 mL 100 mL
<b>Alternate Names:</b>	Plasminogen activator inhibitor-1		
<b>Description:</b>	Plasminogen Activator Inhibitor 1 (PAI-1) depleted plasma is prepared from frozen Mouse plasma using immobilized Anti-Mouse PAI1 IgG.		
<b>UniProt ID:</b>	P22777		
<b>Gene ID:</b>	18787		
<b>Source:</b>	Mouse Plasma		
<b>Formulation:</b>	Frozen plasma with anticoagulant sodium citrate		
<b>Specific Activity:</b>	Undetectable concentration of PAI1 Antigen.		
<b>Storage &amp; Stability:</b>	Store at -80°C for up to 1 year. Upon initial thaw, prepare working aliquots and store at -80°C. <b>Avoid repeated freeze-thaw cycles.</b>		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

