

F10, Sheep Anti-Mouse Coagulation Factor X pAb

Catalog No.	CSI19962A	Quantity:	1.0 ml
	CSI19962B		10 ml

Description: Polyclonal antiserum (host sheep) to mouse factor X.

Concentration: Whole Antiserum

Specificity: Mouse

Host: Sheep

Immunogen: Full length native protein (purified) (Mouse).

Isotype: IgG

Formulation: Liquid

Cross-Reactivity: Mouse

Applications: ELISA, WB

Application Notes: The optimal concentration should be determined by the user for each specific application.

Storage & Stability: When stored at -80°C, product is stable for 5 years from date of delivery. **Avoid repeated freeze-thaw cycles.**

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

