

## F10, Sheep Anti-Mouse Coagulation Factor X pAb

<b>Catalog No.</b>	CSI19963A CSI19963B	<b>Quantity:</b>	1.0 mg 10 mg
<b>Description:</b>	Polyclonal antibody (host sheep) to mouse factor X. IgG fraction purified by immobilized Protein G.		
<b>Concentration:</b>	10.0 mg/ml		
<b>Specificity:</b>	Mouse		
<b>Host:</b>	Sheep		
<b>Immunogen:</b>	Full length native protein (purified) (Mouse).		
<b>Isotype:</b>	IgG		
<b>Formulation:</b>	Liquid in 0.05M Sodium Phosphate + 0.1M NaCl + 1mM EDTA, pH 6.6		
<b>Purification:</b>	IgG fraction purified by immobilized Protein G.		
<b>Cross-Reactivity:</b>	Mouse		
<b>Applications:</b>	ELISA, WB		
<b>Application Notes:</b>	The optimal concentration should be determined by the user for each specific application.		
<b>Storage &amp; Stability:</b>	When stored at -80°C, product is stable for 5 years from date of delivery. <b>Avoid repeated freeze-thaw cycles.</b>		

**NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.**