

## Mouse Anti-Bovine Serum Albumin Clone 12C8D4 mAb

<b>Catalog No.</b>	CSI20161A CSI20161B	<b>Quantity:</b>	1 mg 10 mg
<b>Description:</b>	Mouse Anti-Bovine Serum Albumin Clone 12C8D4 mAb		
<b>Concentration:</b>	10.0 mg/ml		
<b>Host:</b>	mouse		
<b>Immunogen:</b>	Full length native bovine serum albumin (purified)		
<b>Isotype:</b>	IgG1 kappa		
<b>Clone:</b>	12C8-D4		
<b>Formulation:</b>	Frozen liquid buffered in 0.05 M Sodium Phosphate + 0.1 M NaCl + 1 mM EDTA, pH 6.6.		
<b>Purification:</b>	IgG fraction purified by immobilized Protein A.		
<b>Cross-Reactivity:</b>	Cow		
<b>Applications:</b>	ELISA, Western Blot		
<b>Storage &amp; Stability:</b>	Store at -80°C.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.