

Human IgG-Depleted Serum

Catalog No.	CSI20186A	Quantity:	100 mL
	CSI20186B		500 mL

Description: Human serum, depleted of IgG.

Concentration: 0.37 ug/mL Human IgG

Source: Human Serum

Formulation: Frozen Liquid

Purification: Prepared from frozen human serum using immobilized Protein A.

Storage & Stability: store at -80 °C for up to 1 year. **Avoid repeated freeze-thaw cycles.**

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



Cell Sciences®

65 Parker Street

Unit 11

Newburyport, MA 01950

Toll Free: 888-769-1246

Phone: 978-572-1070

Fax: 978-992-0298

E-mail: info@cellsciences.com

Website: www.cellsciences.com