

IgG, Human Plasma, Low Endotoxin

Catalog No.	CSI20191A CSI20191B	Quantity:	50 mg 1 g
Description:	Purified from normal serum by immobilized Protein A using low endotoxin methodology. The purified protein is sterile filtered and lyophilized from azide free PBS. >95 percent pure by SDS-PAGE and <1 EU/mg.		
Concentration:	Lyophilized		
Source:	Human plasma		
Molecular Weight:	160,000		
Formulation:	Lyophilized powder in 0.015M K Phosphate + 0.027M Na Phosphate + 0.15M NaCl, pH 7.2		
Purity:	>95 percent pure by SDS-PAGE and <1 EU/mg.		
Application Notes:	The optimal concentration should be determined by the user for each specific application.		
Storage & Stability:	When stored at -80°C, product is stable for 5 years from date of delivery.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

