

## C57BL6 Mouse IgG, Protein A Purified pAb

<b>Catalog No.</b>	CSI20202A CSI20202B	<b>Quantity:</b>	10 mg 50 mg
<b>Description:</b>	Purified from normal serum by immobilized Protein A. >95 percent pure by SDS-PAGE and preservative free.		
<b>Concentration:</b>	10.0 mg/ml		
<b>Purity:</b>	IgG fraction		
<b>Specificity:</b>	This antibody has been selected as an isotype control as it has no known specificity.		
<b>Host:</b>	C57BL6 Mouse		
<b>Molecular Weight:</b>	160 kDa		
<b>Isotype:</b>	IgG		
<b>Source:</b>	C57BL6 Mouse serum		
<b>Volume:</b>	1.0 ml		
<b>Formulation:</b>	Liquid in 0.02M Sodium Phosphate + 0.15M NaCl, pH 7.4		
<b>Applications:</b>	ELISA, WB		
<b>Application Notes:</b>	The optimal concentration should be determined by the user for each specific application.		
<b>Storage &amp; Stability:</b>	When stored at -80°C product is stable for 5 years from date of delivery. <b>Avoid repeated freeze-thaw cycles.</b>		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.