

CD1 Mouse IgG, Protein A Purified pAb

Catalog No.	CSI20203A CSI20203B	Quantity:	10 mg 50 mg
Description:	Purified from normal serum by immobilized Protein A. >95 percent pure by SDS-PAGE and preservative free.		
Concentration:	10.0 mg/ml		
Purity:	IgG fraction		
Specificity:	This antibody has been selected as an isotype control as it has no known specificity.		
Host:	CD-1 Mouse		
Molecular Weight:	160 kDa		
Isotype:	IgG		
Source:	CD1 Mouse serum		
Volume:	1.0 ml		
Formulation:	Liquid in 0.02M Sodium Phosphate + 0.15M NaCl, pH 7.4		
Applications:	ELISA, WB		
Application Notes:	The optimal concentration should be determined by the user for each specific application.		
Storage & Stability:	When stored at -80°C product is stable for 5 years from date of delivery. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.