

Rabbit IgG, Protein A Purified, pAb

Catalog No.	CSI20228A CSI20228B	Quantity:	500 mg 1 g
Description:	Purified from normal serum by immobilized Protein A. >95 percent pure by SDS-PAGE and preservative free.		
Concentration:	10.0 mg/ml		
Molecular Weight:	160 kDa		
Purity:	IgG fraction		
Specificity:	This antibody has been selected as an isotype control as it has no known specificity.		
Host:	Rabbit		
Isotype:	IgG		
Volume:	50 ml		
Formulation:	Liquid in 0.015M Sodium Phosphate + 0.15M NaCl, pH 7.4		
Source:	Rabbit plasma		
Applications:	ELISA, WB		
Application Notes:	The optimal concentration should be determined by the user for each specific application.		
Storage & Stability:	When stored at -80°C product is stable for 5 years from date of delivery. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.