

## Native Rhesus Monkey IgG

<b>Catalog No.</b>	CSI20237A CSI20237B	<b>Quantity:</b>	10 mg 50 mg
<b>Description:</b>	Protein A Purified Native Rhesus Monkey IgG		
<b>Concentration:</b>	10.0 mg/ml		
<b>Purity:</b>	>95% by SDS-PAGE		
<b>Specificity:</b>	This antibody has been selected as an isotype control as it has no known specificity.		
<b>Host:</b>	Rhesus monkey		
<b>Isotype:</b>	IgG		
<b>Formulation:</b>	Frozen liquid in 0.02 M Sodium Phosphate Buffer + 0.15 M NaCl, pH 7.4. No preservatives and carrier-free.		
<b>Purification:</b>	Rhesus Monkey IgG is purified from normal serum by immobilized Protein A.		
<b>Applications:</b>	ELISA, Western blot		
<b>Storage &amp; Stability:</b>	Store at -80°C. Stable for 5 years from delivery. <b>Avoid repeated freeze-thaw cycles.</b>		

**NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.**