

Sheep IgG Protein G Purified pAb

Catalog No.	CSI20262A CSI20262B	Quantity:	50 mg 100 mg
Description:	Purified from normal serum by immobilized Protein G. >95 percent pure by SDS-PAGE and preservative free.		
Concentration:	4.7 mg/ml		
Host:	Sheep		
Isotype:	IgG		
Conjugate:	Protein G		
Formulation:	Liquid in 0.02M Sodium Phosphate + 0.15M NaCl, pH 7.4		
Purification:	SDS-PAGE and preservative free.		
Applications:	ELISA, WB		
Application Notes:	The optimal concentration should be determined by the user for each specific application.		
Storage & Stability:	When stored at -80°C, product is stable for 5 years from date of delivery. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

