

Polyclonal anti-human Kappa

Catalog No.	CSI20395A	Quantity:	1 mg
	CSI20395B		3 mg

Description: Polyclonal anti-human Kappa light chain procedure: Chickens were immunized with highly purified Human kappa light chains and Fab kappa. After multiple immunizations in Freund's adjuvant eggs were collected, an IgY fraction was prepared and the anti-kappa antibodies were immunoaffinity isolated
This antibody reacts with all human immunoglobulin classes possessing kappa chains and free kappa light chains by IEP and gel diffusion techniques. It does not react with lambda chains.

Concentration: > 1mg/ml

Source: Chicken

Formulation: Liquid, PBS pH 7.2

Purity: Affinity Purified

Storage & Stability: Store at 2-4°C.

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

