

Goat Anti-Human Albumin Affinity Purified pAb

Catalog No.	CSI20403A CSI20403B	Quantity:	1 mg 3 mg
Concentration:	Typically 1mg/ml or greater		
Purity:	Affinity Purified		
Specificity:	This antibody reacts with purified albumin and only with the albumin of human serum.		
Host:	Goat		
Immunogen:	Human Albumin		
Formulation:	Liquid in PBS buffer, pH 7.3		
Applications:	Polyclonal Antibody albumin antibody Working Dilution: Western Blots, 1 to 1,000-10,000 and ELISA, coating, 1 to 200-500		
Storage & Stability:	2-8°C		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

