

Goat anti-Rat CRP Affinity Purified pAb

Catalog No.	CSI20405A CSI20405B	Quantity:	1.0 mg 3.0 mg
Description:	Polyclonal Antibody goat anti-human C-Reactive Protein (anti-CRP) procedure: Goat(s) were immunized with Human CRP C-Reactive Protein . After multiple immunizations in Freund's adjuvant, the goat(s) were bled and serum collected		
Concentration:	Typically 1mg/ml		
Purity:	Affinity Purified		
Host:	Goat		
Immunogen:	Human C-Reactive Protein		
Formulation:	Liquid in PBS, pH 7.2		
Cross-Reactivity:	The affinity purified Polyclonal antibody anti-Human CRP C Reactive Protein reacts only with C-Reactive Protein CRP by ELISA, IEP and Double Diffusion.		
Applications:	Western Blots 1:1,000-1:5,000 ELISA 1:5,000-20,000		
Storage & Stability:	2-8°C		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

