

## Rabbit Anti-Human Ferritin Affinity Purified pAb

<b>Catalog No.</b>	CSI20420A CSI20420B	<b>Quantity:</b>	1 mg 2 mg
<b>Description:</b>	Polyclonal Antibody rabbit anti-human FERRITIN procedure: Rabbit(s) were immunized with purified human ferritin in Freund's adjuvant. The rabbit(s) were bled, serum collected and the anti-ferritin was immunoaffinity purified on an antigen containing immunoabsorbant		
<b>Concentration:</b>	Typically greater than 1 mg/ml		
<b>Purity:</b>	Affinity Purified		
<b>Host:</b>	Rabbit		
<b>Immunogen:</b>	Human Liver Ferritin		
<b>Formulation:</b>	Liquid in PBS, pH 7.2		
<b>Reconstitution:</b>	Polyclonal Antibody rabbit anti-human ferritin Recommended working dilution for: capture antibody: 1:100-1:500. probing antibody: 1:1,000-1:20,000.		
<b>Applications:</b>	IEP and Elisa, the antibodies react with human ferritin.		
<b>Storage &amp; Stability:</b>	2-4°C		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

