

Sheep Anti-Human Elastase pAb

Catalog No.	CSI20424A CSI20424B	Quantity:	0.5 ml 1 ml
Description:	Polyclonal anti-human Elastase procedure: Sheep were immunized with highly purified human elastase. The resulting anti-human elastase was solid phase absorbed to ensure specificity. The sheep Polyclonal anti-human elastase reacts only with human elastase by IEP and gel diffusion techniques		
Host:	Sheep		
Immunogen:	Human Elastase		
Formulation:	Liquid		
References:	Human Neutrophil elastase research: Neutrophil elastase and cancer.Human Neutrophil elastase ref:Surgical Oncology, 22 February 2007 Human Neutrophil elastase research: Association of plasma neutrophil elastase levels with other inflammatory mediators and clinical features in adult patients with moderate and severe pneumonia.Human Neutrophil elastase ref: Respiratory Medicine, 12 February 2007/elastase Human Neutrophil elastase research: Myocardial protection and neutrophil elastase inhibitor: Additive effect of the sivelestat on reperfusion injury.Human Neutrophil elastase ref: Journal of Molecular and Cellular Cardiology, In Press, Corrected Proof, Available online 8 May 2007/elastase		
Storage & Stability:	2-8°C		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.