

Sheep Anti-Human Elastase IgG Fraction pAb

Catalog No.	CSI20425A CSI20425B	Quantity:	1.0 mg 5 mg
Description:	Polyclonal anti-human elastase Procedure: Sheep were immunized with highly purified human elastase. The resulting anti-human elastase was solid phase absorbed to ensure specificity and an IgG fraction was prepared		
Concentration:	Typically > 5mg/ml		
Host:	Sheep		
Formulation:	Liquid in PBS buffer		
Cross-Reactivity:	Polyclonal anti-human elastase Specification: The sheep anti-human elastase reacts only with human elastase by IEP and gel diffusion techniques.		
References:	<p>Human Neutrophil elastase research: Neutrophil elastase and cancer.Human Neutrophil elastase ref:Surgical Oncology, 22 February 2007</p> <p>Human Neutrophil elastase research: Association of plasma neutrophil elastase levels with other inflammatory mediators and clinical features in adult patients with moderate and severe pneumonia.Human Neutrophil elastase ref: Respiratory Medicine, 12 February 2007/elastase</p> <p>Human Neutrophil elastase research: Myocardial protection and neutrophil elastase inhibitor: Additive effect of the sivelestat on reperfusion injury.Human Neutrophil elastase ref: Journal of Molecular and Cellular Cardiology, In Press, Corrected Proof, Available online 8 May 2007/elastase</p>		
Storage & Stability:	2-8°C		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.