

Human Serum Plasma Derived, Pooled Donors

Catalog No.	CSI20428	Quantity:	1 liter
Description:	Plasma derived serum is produced by defibrinating pooled human plasma.		
Source:	Pooled human donors		
Appearance:	Yellow to orange liquid		
Preservatives:	None, sterile filtered, 0.2 micron.		
Country of Origin:	USA		
Storage & Stability:	Store at -20°C to -80°C for 1 year. Upon initial thawing, aliquot and store at -20°C to -80°C. Avoid repeated freeze/thaw cycles.		
Certification:	Non-reactive for HIV-1, HCV, HBV by NAT, HBsAg, HCV Ab, HIV 1 and 2 Ab, and RPR by currently approved FDA methods. However, because no test method can offer complete assurance that infectious agents are absent, this material should be handled at the Bio-safety level 2 (BSL 2) as recommended for any potentially infectious human serum or blood specimen in the CDC/NIH manual "Bio-safety in Microbiological and Biomedical Laboratories", 1999.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



Cell Sciences®
65 Parker Street
Unit 11
Newburyport, MA 01950

Toll Free: 888-769-1246
Phone: 978-572-1070
Fax: 978-992-0298

E-mail: info@cellsciences.com
Website: www.cellsciences.com