

## Goat Anti-Mouse IgG Heavy and Light Chain HRP pAb

<b>Catalog No.</b>	CSI20516A CSI20516B	<b>Quantity:</b>	10 µl 50 µl
<b>Description:</b>	Goat Anti-Mouse IgG Heavy and Light Chain HRP pAb		
<b>Concentration:</b>	1.0 mg/ml		
<b>Specificity:</b>	Mouse IgG Heavy and Light chains		
<b>Host:</b>	Goat		
<b>Immunogen:</b>	Full length native mouse purified protein.		
<b>Isotype:</b>	IgG		
<b>Conjugate:</b>	HRP		
<b>Formulation:</b>	Liquid in 0.01M PBS + 0.0005M EDTA +10% Mannose + 2% Sucrose + 1% BSA + 0.02% Thimerosal, pH 7.6.		
<b>Purification:</b>	Affinity purified and adsorbed against bovine, horse, and human serum proteins.		
<b>Cross-Reactivity:</b>	Normal bovine, horse, and human serum is <2% cross-reactivity by direct solid phase immunoassay.		
<b>Applications:</b>	ELISA		
<b>Application Notes:</b>	ELISA: starting working dilution is 1:10,000. The optimal concentration should be determined by the user for each specific application.		
<b>Storage &amp; Stability:</b>	When stored at 2-4°C, product is stable for 1 year from date of delivery.		

**NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.**