

## Anti-Rabbit IgG, HRP Labeled pAb

<b>Catalog No.</b>	CSI20517A	<b>Quantity:</b>	10 µl
	CSI20517B		50 µl

**Description:** Horseradish peroxidase conjugated goat anti-rabbit IgG, Heavy & Light chain specific. Affinity purified and adsorbed against bovine, horse, human, and mouse serum proteins. Monospecific for rabbit IgG, heavy and light chains, as determined by immunoelectrophoresis against normal rabbit serum. Cross-reactivity with normal bovine, horse, human, and mouse serum is <2 percent by direct solid phase immunoassay. Recommended starting dilution for ELISA is 1:10,000 but proper working dilution must be optimized by the end user

**Concentration:** 160 kDa

**Specificity:** Recognizes Rabbit IgG, Heavy and Light Chains

**Host:** Goat

**Immunogen:** Full length native protein (purified) (Rabbit).

**Purification:** Affinity purified and adsorbed against bovine, horse, human, and mouse serum proteins.

**Cross-Reactivity:** Rabbit

**Applications:** ELISA

**Application Notes:** ELISA: starting working dilution is 1:10,000.  
The optimal concentration should be determined by the user for each specific application.

**Storage & Stability:** When stored at 2-4°C, product is stable for 1 year from date of delivery.

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



**Cell Sciences®**  
480 Neponset Street  
Bldg 12A  
Canton, MA 02021

Toll Free: 888-769-1246  
Phone: 781-828-0610  
Fax: 781-828-0542

E-mail: [info@cellsciences.com](mailto:info@cellsciences.com)  
Website: [www.cellsciences.com](http://www.cellsciences.com)