

## FGA, FGB

### Human Fibrinogen Antigen Capture Plate

<b>Catalog No.</b>	CSI20520A CSI20520B	<b>Quantity:</b>	1 plate 5 plates
<b>Alternate Names:</b>	FGA, FGB		
<b>Description:</b>	Human Fibrinogen Antigen Capture Plate is a 96 well plate with 8 removable strips coated with affinity purified rabbit polyclonal antibody to Human Fibrinogen at 100 ul/well and blocked at 300 ul/well. Intended Use: Specific binding of Human Fibrinogen antigen for sandwich ELISA experiments.		
<b>Samples:</b>	Serum, Plasma, Cell culture supernatants		
<b>Sample Size:</b>	100 µl		
<b>Reagents Included:</b>	1 plate. Other reagents not included.		
<b>Storage &amp; Stability:</b>	Store at 2-4°C. Stable for 1 year from delivery.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

