

## SERPINE1

### Mouse Serpin E1/PAI-1 Antigen Capture Plate

<b>Catalog No.</b>	CSI20532A	<b>Quantity:</b>	1 plate
	CSI20532B		5 plates

**Alternate Names:** PAI1, PAI-1, Planh1, Serpine1

**Description:** Mouse Serpin E1/PAI-1 Antigen Capture Plate is a 96 well plate with 8 removable strips coated with monoclonal antibody to Mouse Serpin E1/PAI-1 at 100 ul/well and blocked at 300 ul/well.

Intended Use: Specific binding of Mouse Serpin E1/PAI-1 antigen for sandwich ELISA experiments.

**Samples:** Serum, Plasma, Cell culture supernatants

**Sample Size:** 100 µl

**Reagents Included:** 1 plate. Other reagents not included.

**Storage & Stability:** Store at 2-4°C. Stable for 1 year from delivery.

**NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.**

