

Plasminogen Depleted Mouse Serum

Catalog No.	CSI20533A CSI20533B	Quantity:	10 ml 50 ml
Description:	Prepared from frozen mouse serum.		
Source:	Mouse Serum		
Purity:	Prepared from frozen mouse serum by immobilized lysine chromatography.		
Storage & Stability:	When stored at -80°C product is stable for 3 years from date of delivery. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



Cell Sciences®
65 Parker Street
Unit 11
Newburyport, MA 01950

Toll Free: 888-769-1246
Phone: 978-572-1070
Fax: 978-992-0298

E-mail: info@cellsciences.com
Website: www.cellsciences.com