

VTN

Mouse Vitronectin Antigen Capture Plate

| | | | |
|---------------------------------|--|------------------|---------------------|
| Catalog No. | CSI20536A CSI20536B | Quantity: | 1 plate 5 plates |
| Alternate Names: | Vn, Vtn | | |
| Description: | Mouse Vitronectin Antigen Capture Plate is a 96 well plate with 8 removable strips coated with affinity purified purified rabbit antibody to Mouse Vitronectin at 100 ul/well and blocked at 300 ul/well. Intended Use: Specific binding of Mouse Vitronectin antigen for sandwich ELISA experiments. | | |
| Samples: | Serum, Plasma, Cell culture supernatants | | |
| Sample Size: | 100 µl | | |
| Reagents Included: | 1 plate. Other reagents not included. | | |
| Storage & Stability: | Store at 2-4°C. Stable for 1 year from delivery. | | |

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

