

Urokinase Coated Plate

Catalog No.	CSI20561A	Quantity:	5 plates
	CSI20561B		25 plates

Description: 96 well plate with 8 removable strips coated with recombinant rabbit Urokinase (uPA) at 100 ul/well and blocked with 300 ul/well. Ideal for irreversible binding of active Plasminogen Activator Inhibitor (PAI-1) for functional ELISA experiments. Rabbit uPA forms a covalent complex with PAI-1 of all species including human, mouse, rat, porcine and rabbit.

Source: Insect cell culture.

Concentration: 0.001 mg/ml

Storage & Stability: Recommended storage 2-4°C. Shelf life 12 months from delivery.

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.