

## Actin, clone 647 Monoclonal Antibody

**Catalog No.:** MON4002

**Quantity:** 1 ml

**Specificity**

This antibody reacts in immunoblots with one single protein band which is corresponding to actin. Antigen origin: cell preparation of chicken skeletal muscle. Antigen location: cytoplasm.

**Immunoglobulin type**

Murine IgM

**Use**

This antibody is useful in immunohistochemistry on frozen sections.  
Not useful for paraffin embedded tissue.

**Instructions for use**

Working dilution of 1:5-1:10 is advised but the optimal dilution should be tested by serial dilution. Positive control: muscle, chicken myofibrils.

**Presentation**

1 ml tissue culture supernatant with 0.1% sodium azide.

**Literature**

- Kaehn, K., et al., 1985, Cell Tissue Res. 239, 417-422.

**FOR RESEARCH USE ONLY, NOT FOR DRUG, DIAGNOSTIC OR OTHER USE.**



**Cell Sciences, Inc.**  
480 Neponset Street  
Bldg 12A  
Canton, MA 02021

Toll Free: 888-769-1246  
Phone: 781-828-0610  
Fax: 781-828-0542

E-mail: [info@cellsciences.com](mailto:info@cellsciences.com)  
Web Site: [www.cellsciences.com](http://www.cellsciences.com)