

Pre-immune Normal Rabbit IgG

Catalog No. NRI01 **Quantity:** 100 µg

Description: This immunoglobulin G (IgG) fraction was purified from normal rabbit serum obtained from animals with no history of immunisation. It can be used as a negative control for purified rabbit IgG antibodies.

Immunogen: none

Isotype: IgG

Formulation: Lyophilized

Purification: Affinity purified from normal rabbit serum.

Reconstitution: **Centrifuge vial prior to opening.** Dissolve in 1 ml sterile water, followed by a buffer suitable to the application..

Storage & Stability: Store as supplied for at least 2 years at 2-8°C. After reconstitution store in working aliquots at -20°C or -80°C. **Avoid freeze-thaw cycles.**

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.